



1kb ladder

Concentration: 104ng/μl

Description

1kb ladder is ideal for determining the size of double-stranded DNA from 500 to 10,000 base pairs. The ladder consists of 10 linear double-stranded fragments. The **2,000bp** and **5,000bp** fragments are present at increased intensity to allow easy identification. All fragments are precisely quantified and mixed during the production. For 5 ul loading, all fragments except 2,000bp and 5,000bp are 40ng. The 2,000bp and 5,000bp fragments are 100ng. This ladder is pre-mixed with loading dye and is ready to use.

Recommended Loading

2-5 μl/Lane

Concentration

Typical Bands 100 ng/5μl

Other Bands 40 ng/5μl

Recommended Electrophoresis Condition

8 cm 1% Agarose Gel, 1×TAE, 7 V/cm, 45min.

Contents (bp)

500、1,000、1,500、**2,000**、3,000、4,000、**5,000**、6,000、8,000、10,000

Storage

1 year at room temperature, for long-time preservation, please storage at -20°C.

